



ENMED® ENM - 212 Disposable Visitor Gown

The ENMED ENM - 212 Visitor Gown is designed to protect from low hazard liquid spray. Visitor Gowns are designed to protect health professionals and patients from cross infection. These gowns have a barrier protection to fluids with the benefit of increased comfort.

Key Features

- · Low lint.
- · Latex free.
- · Fluid repellent.
- · Ultrasonically welded sleeve seams.
- Knitted cuffs for comfort (Optional)

Comfort and Protection



CAT III

Type PB [6] (EN 13688:2013)



Liquid **Protection** Type PB [6] - Light Liquid Splashes (EN 13034:2005+A1:2009)



Biohazard

Type PB [6] - B Infective Agents (EN 14126:2003/AC:2004)

Materials

Fabric: Spunbound Non-woven Fabric

Fabric Weight: 20 GSM Cuff: Knitted cuff (Optional)

Specifications: Anti-Static, Flame Reterdant,

UV Stabilizer

Sewing: Lockstitch sewing for sleeves and shoulders.

Overlock sewing for knitted cuff. Lockstitch for ribbon on the neck.

Belt: 60 cm length, 3 cm width 2 pieces with

Velcro: 5 cm velcro on the neck

Color: Blue/White/Green

Storage and Disposal

- Store in dry, clean conditions in original packaging, away from direct sunlight, sources of high temperature, and solvent vapours
- Store within the temperature range -5°C to +25°C and with relative humidity below 80%
- Expected shelf life is three years from date of manufacture when stored as
- Replace garments if damaged, heavily contaminated or in accordance with local work practice or regulations

Limited Use



DO NOT RE USE

DO NOT DRY

DO NOT WASH

LAUNDERING IMPACTS UPON PROTECTIVE



DO NOT BLEACH

FLAMMABLE



DO NOT IRON



CONSULT FORUSE





CAUTION

Warnings and Limitations

Before use read and understand all user instructions and be sure that the product is suitable for the application and fitted correctly. Product must never be altered or modified. Prior to use inspect the clothing for any damage that could affect its protective function.

Replace any damaged clothing. Care should be taken when removing contaminated garments, so as not contaminate the user with any hazardous subtances. If garments are contaminated, then decontamination procedures should be followed prior to removal of the garment.

The wearing of chemical protective clothing may cause head stress if appropriate consideration is not given to the workplace enviroment. Appropriate under garments should be considered to minimize head stress or damage your garment.





Phone Number +90 324 357 56 44 www.enmedglobal.com

Do not use for:


Contact with heavy oils, sparks or flame, or combustible liquids; Environments with high mechanical risks (abrasions, tears, cuts); Environments with conditions of excessive heat.


Approvals

CE approved under PPE Directive (EU) 2016/425, Category III Certification: Universal Sertifikasyon ve Gözetim Hiz. Tic. Ltd. Şti. Notified Body Numder: 2163

Applications and Performance

Non-Hazardous Particulates	No	Hazardous Dusts & Fibres	No
Non-Hazardous Liquid Splash	No	Hazardous Liquid Splash	No
Non-Hazardous Liquid Spray	Yes	Hazardous Liquid Spray	Yes, if chemical is compatible with suit material
Liquid Continuous Contact/ Immersion	No	Acids/Alkalis	No
Gases and Vapours	No	Organic Solvents	No
Sizing An appropriate size garment should be selected to allow sufficient			





Height

Important Notice

movement for the task.

This guide is only an outline. It should not be used as the only means for selecting protective clothing. Before using any protective clothing, the wearer must read and understand the user instructions for each product. Specific country legislation must be observed. If in doubt, contact a safety professional. Selection of the most appropriate PPE will depend on the particular situation and should only be made by a competent person knowledgeable of the actual working conditions and the limitations of PPE.

Final determination as to the suitability of these products for a particular situation is the employer's responsibility. This information is subject to revision at any time. Always read and follow all User Instructions supplied with your Disposable Protective Gown in order to ensure correct operation. If you have questions contact ENMED Technical Service.

LIMITATION OF LIABILITY: Except as provided above, ENMED shall not be liable or responsible for any loss or damage, whether direct, indirect, incidental, special or consequential arising out of the sale, use or misuse of ENMED products, or the user's inability to use such products.